

EN Instruction manual



FOLD AWAY CHANGING TABLE


AMBER

7734.00


SUMMARY

CONTENTS	PAGES
GENERAL DRAWING	2
-----	-----
DIMENSIONS	2
-----	-----
MODEL	3
-----	-----
TECHNICAL SPECIFICATIONS	3
-----	-----
OPERATION	3
-----	-----
RANGE OF APPLICATION	4
-----	-----
SERVICING	4
-----	-----
DIRECTIONS FOR USE	
1) SIDE RAIL	5
2) ADJUSTABLE BACKREST	5
-----	-----
PERIODIC INSPECTION	6
-----	-----
MAINTENANCE FORM (LOGBOOK)	7
-----	-----

GENERAL DRAWING



DIMENSIONS



MODEL

7734.00 Length 190cm with adjustable backrest and a side rail

TECHNICAL SPECIFICATIONS

FRAME:

- Steel frame, Rilsan coated.

PLATE:

- Support with leatherette cover.

PERFORMANCE:

- Height of the changing table: 80 cm
- Maximum load capacity: 135 kg

WEIGHT:

- 27 kg

OPERATION

The fold away changing table enables the changing of a patient after the shower.

This changing table can be folded away in order to gain room.

The plate remain stable and horizontal thanks to 2 fold-away legs.

The removable side rail secures the patient when on the changing table

The comfort of the patient is improved thanks to the adjustment of the backrest.

RANGE OF APPLICATION

Composed of 2 fold-away legs and of a support with blue leatherette, this changing table is designed to the dressing of a person.

The use of the changing table is reserved to authorized and qualified persons, having carefully read this user guide and aware of the following warnings::

The changing table should be used only for the dressing of a patient.

The maximum weight capacity of the changing table doesn't exceed 135kg

The manufacturer will not be responsible if the devic is used, modified or assembled in another way than prescribed in this user guide!

SERVICING

Do not use abrasive products like « scouring or cleaning powders » to clean the table cover. Use only liquid products (liquid soap).


REVAL advocates the use of BFR 4.0 disinfectant solution to disinfect the changing table (leaflet available on demand)

Please make sure that the table cover is immediately replaced when damaged..

DIRECTIONS FOR USE

1- SECURITY SIDE RAIL


Before each use, please make sure that the 2 pins are in place on the side rail (they are at each end of the side rail).



1- Remove both pins to remove the side rail.

2- Make the side rail revolve upwards. After having put the side rail in place, place the 2 pins in place to block the security rail.

2- ADJUSTABLE BACKREST :



To adjust the slope of the backrest :

- Unscrew the tightening thumb wheel.
- Adjust to the required slope.
- Screw again the tightening thumb wheel (make sure the thumb wheel is correctly tightened to ensure the backrest doesn't fall backwards).

PERIODIC INSPECTION

1) A periodic inspection of the changing table must be done at least once per year. It consists in a visual inspection (in particular of the structure supporting the charge on the table (i.e the fold away legs), as well as security attachment devices (pins on the security side rail) and the thumb wheel tightening the adjustable backrest in the desired position).

2) The periodic inspection has to be completed by someone with the necessary qualifications, having good knowledge of the product and of its use and maintenance. The periodic inspection has to be detailed in the user guide.

3) All important observation has to be notified in the maintenance book which has to be updated by the person in charge of the maintenance / service of the changing table.

The date of corrections brought in answer to the observations has also to be notified in the maintenance book.

4) The person in charge of the inspection will have to notify the date of inspection as well as its result in the maintenance book and sign it.

Elements helping the support of the patient which are detachable will have to be identified and detailed in the maintenance book and controlled. Notifications should bear the conditions of use (nursing home, etc...)

5) In case the periodic inspection would identify defaults, damaged parts or any type of deterioration which could challenge the security of the patient on the fold away changing table, then the owner would have to be informed immediately. If it concerns an immediate danger, then the changing table has to be withdrawn from use. It will be used again only after repair.

6) It is important to notify in the maintenance book any problem or damage occurred between 2 periodic inspections and having already been corrected, challenging the security of the changing table .

7) REVAL SA should be informed of important defaults and damages in order that they take the appropriate action. Forwarded information have to be quoted in the maintenance book.

MAINTENANCE FORM (LOGBOOK)

Control date	Reference of the sling	Tests (135kg)	Observations	Corrections	Agent (name and signature)



FRANCE :

Reval®

Z.A. Fief des Prises

Usseau 17220 Sainte Soulle

Tel: (33) 05.46.42.04.16

Fax: (33) 05.46.43.06.38

www.groupe-reval.com

UNITED KINGDOM :

Reval® Continuing Care Limited

Unit 17 Brynmenyn Busines Centre

St Theodore's Way - Brynmenyn

CF32 9TZ BRIDGEND - WALES

Tel: 01656 721 001

Fax: 01656 720 006

www.reval-group.com