


Standing Harness

The Oxford® Standing Harness is a versatile sling to assist with standing, walking or toileting. It is not a general-purpose sling and requires the patient to be able to partially weight-bear. Often used for rehabilitation activity in combination with an overhead lift.



- ✓ Colour coded binding indicates sling size
- ✓ Coloured attachment loops promote correct/safe positioning
- ✓ Padded body and harness straps for comfort
- ✓ Fully adjustable belt and straps to meet patient needs
- ✓ Can be customised to fulfil patient specific needs




Chair to Chair



Bed to Chair/Chair to Bed
(From/to seated position)



Chair to Commode/Toilet



Transfer to Bath



Floor to Bed/Chair



Rehab/Assisted Standing

Size	Part No.	SWL
XS	SL26600	227 kg
S	SL26601	227 kg
M	SL26602	227 kg
L	SL26603	227 kg
XL	SPECIAL	227 kg



Joerns Healthcare strongly recommends that a thorough risk assessment is conducted by a competent person prior to lifting/transferring any client to ensure their complete safety, comfort and dignity is maintained at all times.