



Electrical all functions psychiatric bed

# Sana 1200

Bed has been designed with the psychiatric hospital caregivers' help



BVCert. 6033000



## Patients' benefits

- **Patient's safety :**

- High resistance sleep deck, equipped with welded slats.
- Resistant structure that will keep its shape under all conditions
- All the moving parts of the bed are strongly fixed. Bed cannot be dismantled by dementia patients.
- Patient immobilisation : Bed provides contention capability; the patient's movements can be limited using safety straps.

- **Emergency care :** The bed is equipped with a damped CPR function.

- **Sleep deck extension :** An 18cm strongly secured frame extension is available.

- **Patient Safety :** For therapeutic needs, easy fastening folding siderails can be installed.



Fixe Resistant structure with slats welded frame



Head section fitted with a gas spring  
•Knee section fitted with a gas spring.



Articulated sleep deck with a head section fitted with a gas spring

Mechanical leg section or fitted with a gas spring

Mechanical tibia section



## Caregivers' benefits

- **Motorisation :** The electrical functions movements ensures the patient and caregiver's safety and comfort.

- **Mancœuvrability :** Due to 125mm individual brake casters, patient transfer is easy.

- **Safety optimization :** Wireless control hand pendant available (as option) to ensure dementia patient safety. In case of need the cable can be removed and the bed can be powered by battery alone.

- **Safety increased :** Bed is equipped with a selective lock out of electric bed functions.

- **Fastenings are available to receive any therapeutic accessories :** patient helper, IV pole, ...

- **Optimal wall protection :** A metal spacer keeps the bed sufficiently away from the wall to prevent damage when the bed is being raised, or in a chair position

- **Cleaning optimization :** Easy to clean and disinfect . All the materials are compatible with standard chemical cleaning products.





Electrical all functions psychiatric bed

# Sana 1200

Bed has been designed with the psychiatric hospital caregivers' help



BVCert. 6033000

## Technical specifications

### Key features

- Frame with standard sleep deck : 200 x 90 cm
- Maximum width with siderails: 102 cm
- Maximum length with siderails : 211 à 219 cm
- Sleep deck height position : 38 à 83 cm
- Head section angle : 70°
- Knee section angle: 30°
- Foot section angle : 20° à 30°
- Safe working load: 185 kg
- Bed lift maximum capacity: 250 kg

### Standard equipment

- Manual head section
- Electrical hi-lo
- Central locking from the bed foot
- hand pendant control
- 2 Patient Helper supports from the head of the bed
- 4 I.V. Pole supports

### Option

- Electrical head section
- Electrical knee section with a manual foot section
- 4 angle bumpers
- Integrated supports for folding siderails
- Wireless hand pendant control
- Possibility to replace casters with floor fixes
- Rotating Patient Helper with plastic handle and adjustable strap
- Fixed or adjustable height I.V. Pole

