



Fixed height psychiatric bed
with assisted movements

Sana 400

Bed has been designed with psychiatric hospital caregivers



BVCert. 6033000

Patients' benefits



- **Patient's dignity kept** : The comfort chair position allows the patient to keep his dignity while interacting with people.
- **Patient safety** :
 - Resistant structure that will keep its shape.
 - All the moving parts of the bed are strongly secured and it cannot be dismantled by dementia patients.
- **Patient safety** : For therapeutic needs, optional easy fastening folding siderails can be installed.
- **Patient immobilization** : the patient's movements can be limited by using safety straps.
- **Sleep deck extension** : An 18cm strongly secured frame extension is available as an option.



Caregivers' benefits

- **Cleaning optimization** : All the materials are easy to clean and resist the use of standard chemical cleaning products.
- **Manoeuvrability** : 125mm Individual brake casters are available on this fixed height bed
- **Accessories** : Fastenings are available to receive any therapeutic accessories: patient helper, IV pole





Fixed height psychiatric bed with assisted movements

Sana 400

Bed has been designed with psychiatric hospital caregivers



BVCert. 6033000

Technical specifications

Key features

- Frame with standard sleep deck : 200 x 90 cm
- Maximum width with siderails : 102 cm
- Maximum length with boards : 211 à 219 cm
- Sleep deck height 38 cm
- Head section angle 70°
- Knee section angle: 30°
- Foot section angle : 20° à 30°
- Safe working load : 185 kg
- Bed lift maximum capacity: 250 kg

Standard equipment

- Fixed height sleep deck
- Central locking from the foot of the bed
- 2 Patient helper supports integrated in the head of the bed

Options

- Manual head section assisted by a gas spring
- Manual Knee section assisted by a gas spring
- Manual foot section
- Bumpers located at 4 corners
- Integrated supports for folding siderails
- Possibility to replace casters by floor fixation
- Rotating Patient helper with a plastic handle and an adjustable strap
- Fixed or adjustable I.V. Pole

