



Hydraulic secured psychiatric bed



Sana 800

Bed has been designed with psychiatric hospital caregivers' help



BVCert. 6033000

Patients' benefits

- **Patient's dignity is always maintained** : due to the comfort chair position.
- **Patient safety** :
 - Resistant structure that will keep its shape under all conditions
 - All the moving parts of the bed are strongly secured: the bed cannot be dismantled by dementia or psychiatric patients.
- **Ergonomy** : The hydraulic hi-low movements ensure comfort and safety of both patients and caregivers.
- **Patient immobilisation capability** : the patient's movements can be limited by using safety straps.
- **Sleep deck extension** : An 18cm strongly secured frame extension is available.
- **Emergency care** : The bed is equipped with a damped CPR function.
- **Patient safety** : For therapeutic needs, optional easy fastening folding siderails can be installed.

Caregivers' benefits

- **Cleaning optimization** : all the materials are easy to clean and resistant to the use of standard chemical cleaning products.
- **Manœuvrability** : Due to 125mm individual brake casters, patient transfers around the hospital are available.
- **Accessories** : Fastenings are available to receive any therapeutic accessories: patient helper, IV pole,
- **Optimal walls protection** : A metal spacer keeps the bed sufficiently away from the wall to prevent damage when the bed is being raised or in the chair position.



Fixe Resistant structure with slats welded frame



Head section fitted with a gas spring
•Knee section fitted with a gas spring.

Mechanical tibia section



Articulated sleep deck with a head section fitted with a gas spring

Mechanical leg section or fitted with a gas spring



Hydraulic secured psychiatric bed

Sana 800

Bed has been designed with psychiatric hospital caregivers' help



BVCert. 6033000

Technical specifications

Key features

- Frame with standard sleep deck : 200 x 90 cm
- Maximum width with siderails : 102 cm
- Maximum length with siderails : 211 à 219 cm
- Sleep deck mini and max height position : 38 à 83 cm
- Head section angle : 70°
- Knee section angle: 30°
- Foot section angle : 20° à 30°
- Working safe load: 185 kg
- Bed lift maximum capacity: 250 kg

Standard Equipment

- Hydraulic hi-lo
- Central locking from the foot of the bed
- 2 supports for patient helper located at the head of the bed
- 4 supports for iv poles

Options

- Manual head section assisted by a gas spring
- Manual knee section assisted by a gas spring
- Manual foot section
- Bumpers located at 4 corners
- Integrated folding siderails
- Wireless hand pendant control
- Possibility to replace casters by floor fixation
- Rotating Patient helper with a plastic handle and an adjustable strap
- Fixed or adjustable I.V. Pole

