


Quickfit

The Oxford® Quickfit sling is an easy to fit, general-purpose sling with divided leg sections, designed to suit around 85% of clients. The Quickfit design has become one of the most recognised sling choices available today.


- ✓ Available in polyester solid or net material options
- ✓ Optional padded leg sections for increased comfort
- ✓ Colour coded binding indicates sling size
- ✓ Positioning handles allow easy patient repositioning
- ✓ Coloured attachment loops promote correct/safe positioning
- ✓ Modesty loop helps protect patient dignity during transfer
- ✓ Can be customised to fulfil patient specific needs




Chair to Chair



Bed to Chair/Chair to Bed



Chair to Commode/Toilet



Transfer to Bath



Floor to Bed/Chair
(Head Support recommended)



Rehab/Assisted Standing

Size	Part No. (Solid)	Part No. (Net)	Part No. (Padded)	SWL
XS	SPECIAL	SPECIAL	SPECIAL	227 kg
S	SL1000	SL1001	SL1000P	227 kg
M	SL1003	SL1004	SL1003P	227 kg
L	SL1006	SL1007	SL1006P	227 kg
XL	SPECIAL	SPECIAL	SPECIAL	227 kg



Joerns Healthcare strongly recommends that a thorough risk assessment is conducted by a competent person prior to lifting/transferring any client to ensure their complete safety, comfort and dignity is maintained at all times.