


Comfort

The Oxford® Comfort sling is an easy fit, contoured sling designed to accommodate the majority of clients. For use in conjunction with a 4-point positioning cradle, the patented Securi3 clip system ensures no inadvertent detachment of the sling.

- ✓ Available in polyester solid or net material options
- ✓ Optional padded leg sections for increased comfort
- ✓ Patented Securi3 clip attachment points
- ✓ Colour coded binding indicates sling size
- ✓ Positioning handles allow easy patient repositioning
- ✓ Modesty loop helps protect patient dignity during transfer
- ✓ Can be customised to fulfil patient specific needs




Chair to Chair



Bed to Chair/Chair to Bed



Chair to Commode/Toilet



Transfer to Bath



Floor to Bed/Chair



Rehab/Assisted Standing

Size	Part No. (Polyester)	Part No. (Net)	Part No. (Padded)	SWL
XS	SL25501	SL25511	SL25506	227 kg
S	SL25502	SL25512	SL25507	227 kg
M	SL25503	SL25513	SL25508	227 kg
L	SL25504	SL25514	SL25509	227 kg
XL	SL25505	SL25515	SL25510	227 kg



Joerns Healthcare strongly recommends that a thorough risk assessment is conducted by a competent person prior to lifting/transferring any client to ensure their complete safety, comfort and dignity is maintained at all times.