


Deluxe Standing Disposable

The Oxford® Disposable Standing sling is an easy fit sling suitable for standing and supporting. It is not a general purpose sling and clients should be assessed as to their suitability. It will not suit double amputees or lifting from a horizontal position. In the interests of avoiding cross-infection, this sling is specifically designed for single patient use only.



- ✓ Fully padded design offers excellent patient comfort
- ✓ Coloured attachment loops promote correct positioning
- ✓ Non-slip back section ensures the sling does not 'ride up'
- ✓ Adjustable waist strap for added support and security




Chair to Chair



Bed to Chair/Chair to Bed
(From/to seated position)



Chair to Commode/Toilet



Transfer to Bath



Floor to Bed/Chair



Rehab/Assisted Standing

Size	Part No.	SWL
XS	SPECIAL	200 kg
S	SL13525N	200 kg
M	SL13526N	200 kg
L	SL13527N	200 kg
XL	SPECIAL	200 kg



Joerns Healthcare strongly recommends that a thorough risk assessment is conducted by a competent person prior to lifting/transferring any client to ensure their complete safety, comfort and dignity is maintained at all times.